GIMA

FONTE DI LUCE 150W
LIGHT SOURCE 150W
SOURCE LUMINEUSE 150W
FUENTE DE LUZ 150W
FONTE DE LUZ 150W
ΠΗΓΗ ΦΩΤΟΣ 150W

MANUALE D'USO E MANUTENZIONE
USE AND MAINTENANCE BOOK
INSTRUCTIONS DE FONCTIONNEMENT ET ENTRETIEN
MANUAL DE USO Y MANTENIMIENTO
MANUAL DE USO E MANUTENÇÃO
ΕΓΧΕΙΡΙΔΙΟ ΧΡΗΣΗΣ ΚΑΙ ΣΥΝΤΗΡΗΣΗΣ

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
CE mark: 93/42 EEC directive

Medical device: class I

SYMBOLS	
BF Type apparatus	
CE mark - 93/42 EEC directive	
Fuses	
Alternating current	

STANDARD COMPONENTS

GIMA 150 Watt Light source

Main cable

Halogen lamp 150 Watt

User manual

FEATURES

Gima Light source is complete with a 150 Watt halogen bulb. Light emission is variable from 50 to 150 W by a luminous intensity regulator. A continuous forced ventilation system produces effective cooling of the unit. The luminous intensity regulator mounted on the unit is a rheostat type and allows energy saving and a longer life for the bulb. The transformer is certified and it has a 100°C thermal protection.

The instrument must be operated by qualified personnel only and must only be used for those applications for which it was specifically designed, any other use is to be considered inappropriate and dangerous and the producer declines any responsibility for damages originated by inappropriate use.

PRESCRIPTIONS

Make sure the unit is grounded to assure correct electrical safety of the instrument in accordance with the present electrical protection standards.

The producer declines any responsibility for any eventual damage originated by an inappropriate grounding of the instrument.

Never try to open, modify or repair the instrument, but contact GIMA technical service.

Lamp or fuses replacement (following instruction at page 2) is possible only by qualified personnel.

Never touch the product with wet hands or feet

Never use the product without wearing shoes

Never leave the product exposed to atmospheric phenomena

Never leave the product on when not in use, always switch off

Never use or turn on the product in presence of flammable gas

Light source packed for transport and storage may be exposed to the following environment at conditions:

- 1) Ambient temperature between -40°C and +70°C
- 2) Relative humidity between 10% and 100%, including condensation
- 3) Atmospheric pressure between 500 and 1060 hPa

FUNCTIONING

- 1) Unpack the product and check everything is in order, if in doubt do not use the unit and call GIMA.
- 2) Plug the unit into a 220-230V 50 Hz. power source (use only the cable supplied with the unit, when using extensions or adaptors make sure these are fitted with the ground conductor).
- 3) Switch on the Light source setting the main switch on the front to „ON“, a light just above the main switch and the 150 watt light will light up.
- 4) Screw the GIMA Light source adaptor to the cable (cable and adaptors are not standard components of the the unit, but to use the light is necessary to have cable and adaptor)
- 5) Connect the cable to GIMA Light source
- 6) Using the knob select the desired light intensity, from 50 to 150 watt, move right to increase intensity, move left to decrease intensity

MAINTENANCE

GIMA Light source does not need any special maintenance. Just clean out dust and dirt with a soft alcohol impregnated cloth. Keep GIMA Light source in a place with normal climate, in following ratings:

- A) ambient temperature between 10°C and 40°C
- B) relative humidity between 30% and 75%
- C) atmospheric pressure between 700 and 1060 hPa.

Avoid any repair or alteration of device by unqualified and untrained persons.

FUSES REPLACEMENT

In case of a blown fuse proceed as follows:

- 1) Remove power cable from the socket
- 2) Open the fuseholder on the back side of the instrument
- 3) Replace the fuses with new ones

REPLACEMENT OF THE LAMP

- 1) Take out power cable from the socket
- 2) Wait for few minutes for the 150 watt light to cool
- 3) Unscrew the closing panel on the light source lid
- 4) Pull up the lamp holder
- 5) Take out the lamp with a light pressure down and upwords GIMA Light source works only with original GIMA lamps (15Volt 150W). Do not touch the mirror of the lamp with your fingers, you can seriously damage it causing power leakage

GUARANTEE

Thank you for having acquired our product.

This product meets the most stringent requirements regarding the selection of quality manufacturing materials and also the final control. The product has a 12-month warranty, valid from the date it is delivered by GIMA.

During the guarantee period free repair and/or replacement of any defective parts due to faulty manufacture will be given, labour, postal, transport and packaging charges etc. are not included. The guarantee therefore excludes components subject to wear and tear such as parts in rubber or PVC, doppler probes, SpO₂ sensors, lamps, batteries, electrodes, handpieces, washers and packings, resistances, etc.

No compensation will be given for loss of use of the product.

Furthermore replacement or repairs effected during the guarantee period do not lengthen the duration of the guarantee.

This guarantee is void in the event of: repairs effected by unauthorised persons or with spare parts not approved by GIMA, damage or defects caused by negligence, blows, abnormal use of the appliance or faulty installation.

The guarantee will be invalidated if the registration number has been removed, cancelled or altered. Faulty appliances must be returned only to the retailer where the item was bought. Any item dispatched directly to us will be refused.

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

Disposal: The product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment. For further information on recycling points contact the local authorities, the local recycling center or the shop where the product was purchased. If the equipment is not disposed of correctly, fines or penalties may be applied in accordance with the national legislation and regulations.

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