CINTURA DI SICUREZZA
SAFETY BELT
SANGLE DE SÉCURITÉ
CINTURÓN DE SEGURIDAD
CINTO DE SEGURANÇA
SICHERUNGSGURT
ΖΩΝΗ ΑΣΦΑΛΕΙΑΣ
PAS BEZPIECZEŃSTWA
BEZPEČNOSTNÍ PÁS
SÄKERHETSBJÄLTA
TURVAVYÖ
VARNOSTNI PAS
BEZPEČNOSTNÝ PÁS
CENTURĂ DE SIGURANȚĂ
VEILIGHEIDSGORDEL
SIGURNOSNI POJAS
BIZTONSÁGI HEVEDER
SIKKERHEDSSELE
ΠΡΟΤΡΑΞΗ ΚΟΛΑΗ
SAUGOS DIRŽAS
DROŠĪBAS JOSTA
TURVAVÖÖ

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REF B2 (GIMA 33840) - B4 (GIMA 33842) - B8 (GIMA 33843) - B9 (GIMA 33844) - B11 (GIMA 33847) - B12 (GIMA 33849)

JIANGSU RIXIN MEDICAL EQUIPMENT CO.,LTD.
No.427 Yangjin Road, Jinfeng, Zhangjiagang, Jiangsu
Province, CHINA
Made in China

SUNGO Europe B.V., Fascinatio Boulevard 522, Unit 1.7,
2909VA Capelle aan den IJssel, The Netherlands

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

ENGLISH

INSTRUCTIONS FOR USE

Name

Safety Belt

Model

B2, B4, B8, B9, B10, B11, B12

Intended Use

Safety Belt is used for fix the patient slips and falls from the stretcher. The Safety Belt is used with the stretcher, provided non-sterile.

Component

Safety Belt is mainly composed of Polyester straps and buckle.

Material

Safety Belt is made of Polypropylene or Polyester

Specification

The dimension is as below:

Model	Specification	Picture	Material
B2	180cm		Polypropylene
B4	180cm		Polyester
B8	180cm		Polyester
B9	180cm		Polyester
B10	110cm		Polyester
B11	180cm		Polypropylene
B12	166cm		Polypropylene

Storage

Store in moisture-proof and non corrosive environment.

Cleaning and maintenance

Keep clean as routine. Wipe with 75% alcohol solution.

Transportation

The Safety Belt can use general means of transport.

No inversion is allowed during storage and transportation, and stacking is carried out as shown in the packing box.

Shelf Life

8 years

Use method

Model	Picture	Use method
B2		To secure: Weave straps through handhold on backboard, and around patient and board, couple the male end into the female end until you hear a click sound, then properly buckle. To release: Press side of closure, then pull apart the two ends.

B4		To secure: Weave straps through handhold on backboard, and around patient and board, couple the male end into the female end until you hear a click sound, then properly buckle. To release: Press side of closure, then pull apart the two ends.
B8		To secure: Put belt's two hooks on backboard pin, and around patient and board, couple the male end into the female end until you hear a click sound, then properly buckle. To release: Press side of closure, then pull apart the two ends.
B9		To secure: Weave straps through handhold on backboard, and around patient and board, couple the ends by pushing metal tongue into buckle slot until you hear a click sound, then properly buckle. To release: press down on the button centrally located on top of the buckle and pull metal tongue out of buckle slot.
B10		1. Wrap webbing around stretcher frame or around handhold in backboard and through end of webbing where loop is sewn. 2. Put the two belts ends between the buckle and metal tongue, couple the ends by pushing metal tongue into buckle slot until you hear a click sound, then properly buckle. 3. To release: press down on the button centrally located on top of the buckle and pull metal tongue out of buckle slot.
B11		To secure: Weave straps through handhold on backboard, and around patient and board, couple the ends by pushing metal tongue into buckle slot until you hear a click sound, then properly buckle. To release: press down on the button centrally located on top of the buckle and pull metal tongue out of buckle slot.
B12		1. Weave straps through handhold on backboard, and around patient and board. 2. Shoulder straps should around patient shoulder and head side on board.

Precaution and Warning

1. Always keep clean (including disinfection).
2. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Indice dei simboli - Symbol index - Index des symboles - Índice de símbolos - Índice de símbolo - Symbolindex - Ευρετήριο συμβόλων - Indeks symboli - Index symbolů - Symbol index - Symboli-indeksi - Indeks symbolov - Index symbolov - Index de simbol - Symbol index - Indeks simbola - Szimbólum index - Symbolindeks - Индекс на символа - Simbolių rodyklė - Simbolu rādītājs - Símbolite indeks

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