PROFESSIONAL MEDICAL PRODUCTS

IMMOBILIZZATORE SPINALE
SPINAL IMMOBILIZER
IMMOBILISATEUR DE COLONNE VERTÉBRALE
INMOVILIZADOR ESPINAL
IMOBILIZADOR ESPINHAL
WIRBELSÄULEN-WEGFAHRSPERRE
СПИНАЛЕН ИМОБИЛАЙЗЕР
IMOBILIZÉR PÁTEŘE
SPINAL IMMOBILISER
ΑΚΙΝΗΤΟΠΟΙΗΣΗ ΣΠΟΝΔΥΛΙΑΣ
SPINAL IMMOBILISER
SPINALNI IMOBILIZATOR
GERINCS INDÍTÁSGÁTLÓ
SPINALINIS IMOBILIZATORIUS
MUGURA IMOBILIZERS
RUGGENGRAALIMMOBILISATOR
IMOBILIZER KRĘGOSŁUPA
IMOBILIZATORUL COLONIALĂ
IMOBILIZÉR CHRBTICE
SPINALNI IMOBILIZATOR
SPINAL IMMOBILISERARE

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REF BS-01 (GIMA 34032) - BS-01 (GIMA 34062)

JIANGSU RIXIN MEDICAL EQUIPMENT CO.,LTD.
No.427 Yangjin Road, Jinfeng, Zhangjiagang, Jiangsu
Province, CHINA
Made in China

SUNGO Europe B.V., Fascinatio Boulevard 522, Unit 1.7,
2909VA Capelle aan den IJssel, The Netherlands

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

ENGLISH

INSTRUCTIONS FOR USE

Name: Splint

Model

PS-01, PM-01, BS-01, BS-02, AS-01, AS-02, AS-03, AS-04, AS-05A, AS-05B, AS-05C, AS-05D, AS-05E, AS-05F, TS-01, TS-02, AM-01

Intended Use

Splint is used for fix each part of the body.

Component

Splint is mainly composed of frame, foam and burrs and zipper, safety belt, air pump and valves.

Material

Splint is made of Aluminum Alloy, PVC, EVA, TPU and Oxford leather.

Specification

The dimension is as below:

Model	Specification	Picture	Material
BS-01 (size M) BS-01 (size L)	82x79cm 88x86cm		PVC

Storage

1. Store in moisture-proof and non corrosive environment.
2. Stacking and stacking according to the boxes in storage and transportation.

Cleaning and maintenance

Keep clean as routine. Wipe with 75% alcohol solution. Check the connecting and rotating parts of the stretcher. Check the safety belts regularly, if they are worn or have cracks, replace them immediately.

Transportation

The Splint can use general means of transport. No inversion is allowed during storage and transportation, and stacking is carried out as shown in the packing box.

Shelf Life

8 years

Use method

Model	Picture	Use method
BS-01 (size M/L)		<p>The device is used for fixed adult.</p> <ol style="list-style-type: none"> 1. Put the splint at the back of the patient when the patient is sitting. Close up the fixing pieces for body and connect them through fixing belts and buckles. Fold and fix the rest belts. If the patient is injured in cervical, first fix the patient with a cervical collar, then wrap up patient's head by the fixing pieces (use together with cushion if needed, put the cushion at the back of the patient before use), fix patient's forehead and chin by fixing belts. 2. After use, clean, disinfect and dry the splint, fold it according to the way before folding, then put into storage.

Precaution and Warning

1. Always keep clean (including disinfection).
2. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

GIMA WARRANTY TERMS

The Gima 12-month standard B2B warranty applies.

Indice dei simboli - Symbol index - Index des symboles - Índice de símbolos- índice de símbolo - Symbolindex - Ευρετήριο συμβόλων - Indeks symboli - Index symbolů - Symbol index - Symboli-indeksi - Indeks symbolow - Index symbolow - Index de simbol - Symbol index - Indeks simbola - Szimbólum index - Symbolindeks - Индекс на симбола - Simbolių rodyklė - Simbolu rādītājs - Sümboleite indeks

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