

iHEALTH WIRELESS PERSONAL HEALTHCARE PRODUCTS


iHEALTH MYVITALS 2.0 APP - free

To manage wireless:

- 1 iHEALTH B.P. MONITOR - ARM
- 2 iHEALTH B.P. MONITOR - WRIST
- 3 iHEALTH SCALE
- 4 iHEALTH BODY ANALYSIS SCALE (9 parameters)
- 5 iHEALTH BODY ACTIVITY TRACKER - WRIST/BELT
- 6 iHEALTH PULSE OXIMETER
- 7 iHEALTH GLUCOSE MONITOR



iHEALTH products to manage your personal health
 iHealth designs innovative, mobile personal healthcare products for everyone. iHealth wireless products with **Bluetooth** connectivity for iPad, iPhone, iPod and Android system, allow to easily measure and track personal data and share instantly with family or doctor. The powerful **iHealth MyVitals 2.0 App** will also allow to manage all your health data with custom graphs, statistics, trends. This unique application **tracks the relation between various parameters** such as blood pressure, weight, calorie intake and activities. App compatible with all Apple and Android platforms. Includes a free iHealth Cloud service. MyVitals App available in: GB, FR, IT, ES, DE, TR, JP, KR, CN. Other EU languages available in next months. Product package and manual in 5 languages (GB, FR, IT, ES, DE). On request GB, NO, SE, FI, DK, DE. All products have 2 years warranty.



1 2 iHEALTH B.P. MONITOR - ARM/WRIST



TECHNICAL SPECIFICATIONS		
	BP5	BP7
Classification internally	Type B applied part, IPX0 no AP or APG	Type BF
Wireless communication:	Bluetooth V3.0+EDR	Bluetooth V3.0+EDR class 2 SPP
Size:	145x58x30 mm	72x74x17.6 mm
Cuff circumference:	22-42 cm	13.5-22
Weight:	approx 135 g (excluding cuff)	105 g
Measuring method:	oscillometric, automatic inflation and measurement	
Power:	DC: 5 V - 1.0A, Battery: 1x3.7 V - Li-ion 400 mAh	
Cuff pressure range:	0-295 mmHg ±3 mmHg	0-300 mmHg ±3 mmHg
Pulse rate range:	40-180 beats/min, accuracy; ±5%	
Battery life:	300 recharges	
Numbers of readings with charged battery:	120 readings	250 readings



- 23500 BP5 WIRELESS ARM BLOOD PRESSURE MONITOR
- 23501 BP7 WIRELESS WRIST BLOOD PRESSURE MONITOR
- measure blood pressure and heart rate directly from tablets or mobile devices
- record systolic/diastolic number, heart rate and pulse wave
- multiple users and unlimited memories
- view historical readings and WHO classifications (only BP5)
- generate custom graph by date or BP numbers

3 iHEALTH SCALE

BMI function in the APP




- 23505 HS3 WIRELESS SCALE
- 23506 HS4 WIRELESS SCALE
- A high quality design digital scale, with slim lightweight design and reliable tempered glass platform, the iHealth scales works online and offline. Thanks to the iHealth MyVitals App it is possible to view single readings or long term trends. Track calories and daily activities, set target goals for weight management, organize records with personalized graphs and unlimited memories.

TECHNICAL SPECIFICATIONS		
	23505	23506
Max weight capacity:	150 kg (330 lbs)	180 kg (400 lbs)
Size:	330x330x39 mm	350x350x28.5
Accuracy:	±1%	±1%
4x1.5V AAA batteries included		



4 iHEALTH BODY ANALYSIS SCALE

TECHNICAL SPECIFICATIONS	
Classification:	Internally powered, Type BF applied part
Wireless Communication:	Bluetooth V3.0+EDR Class 2 SPP
Size:	380x310x35 mm (14.96" x 12.20" x 1.38")
Weight:	1.5 kg (3.3 lbs) (excluding battery)
Measuring method:	Automatic full electronic measurement
Power:	4 x 1.5 V AA batteries
Body weight range:	5-150 kg (11-130 lbs)
Body fat measurement range:	5%-65%
Body water measurement range:	20%-85%
Visceral fat rating range:	1-59
Body weight:	±0.5 kg (5 ~ 40 kg); ±1%+0.1 kg (40 ~ 150 kg)
Body fat:	±1%
Body water:	±1%
Body muscle mass:	±1% + 0.1 kg (±1% + 0.2 lbs)
Bones mass:	±0.3 kg (0.66 lbs)
Visceral fat rating:	±2
DCI (Daily Calorie Intake):	±200 kcal
Switchable kg/lb	



- 23507 HS5 WIRELESS BODY ANALYSIS SCALE - 9 parameters
- Measures 9 key data points: **weight, body fat, lean mass, muscle mass, bone mass, body water, daily calorie intake, body mass index and visceral fat rating.** Log calories and daily activities, plan weight goals. Organize records, charts, historical readings and view trends. Unlimited users and memory (offline memory of 200 measures). Possibility to set measurements and activity reminders.

