

■ INFRARED THERAPY LAMP - 250 W

- 28653 INFRARED THERAPY LAMP - 250 W - trolley
- 28654 SPARE BULB - 250 W

The infrared light effectively relieves the suffering caused by sore muscles and joints ankylosis. The comfortable warmth created by infrared light penetrates deep into the skin, stimulates blood circulation and warms the muscles. While the muscles are relieved, they automatically stretch and relax.

At the same time, the blood circulation helps your body get rid of impurities and quickly sends oxygen-rich blood to stressed or aching muscles, causing an effective pain relief.



■ INFRAREX: THERAPEUTIC HEATING DEVICE

- 28650 INFRAREX 230 V - 50/60 Hz
- 28651 SPARE BULB 7.5 W



Infrarex is a modern scientifically electronic therapeutic heating device which provides deep, infrared heat that can penetrate far below the surface of the skin to produce a shooting, comforting effect upon the nerves. It relieves muscular tension, reduce inflammation and provides pain relief without medication for many conditions such as: arthritis, bursitis, sinusitis, tennis elbow, back ache, neuralgia, muscular pains, aches, sprains, and other musculoskeletal aches.

Infrarex provides pain relief without medication, you will experience pain relief within 10-15 minutes. No concern for side effects, simple to be used.



■ ELECTRIC HEATING PAD

- 28655 ELECTRIC HEATING PAD Suitable for many applications.
- 28656 ELECTRIC "SAND" HEATING PAD

The use of purified marine sands make it possible to naturally transmit their beneficial warmth.

Heating: 48-52°C (initial 30 minutes) 52-56°C first hour, up to 65°C following hours.


External cover in cotton (94%) and nylon, internal bag in PVC
Operating power: 230 V - 50/60 Hz - 75 W (28656) - 55 W (28655)



■ ACCESSORIES FOR ANY BRAND OF NEBULIZER

ACCESSORIES FOR NEBULIZER

Gima offers different nebulizer bulbs with their line of accessories, suitable for any brand of nebulizer.



FASTERJET BULB


with Valve System

A new designed Jet Nebulizer capable of delivering any medication in a very short time with a high amount of respirable particles. The Valve System makes therapy extremely efficient by reducing the amount of wasted medication while optimizing the output to the patient.

TECHNICAL SPECIFICATIONS FASTERJET

Capacity: 16 ml Neb rate: 0.60 ml/min
Respirable volume: <math>< 5 \mu\text{m}</math> 80%
Residual volume: 0.8 ml
MMAD (Mass Median Aerodynamic Diameter): 2.7 μm

Fasterjet bulb:
For delivering any medication in a very short time



GIMA code	NEBULIZER ACCESSORIES	To be used with
Standard accessories with Famiglia, Corsia, Eolo		
28158	PVC connection tube - 1m	all nebulizers
28139	Nebulizer bulb	all nebulizers
28146	Facemask - adult* - replaces 28143	bulb 28139/28140
28147	Facemask - paediatry* - replaces 28144	bulb 28139/28140
28157	Mouthpiece - nasal prong	bulb 28139/28140
28138	Nebulizer kit <i>adult and pediatric facemask, pvc-tube, bulb, nasal prong, mouthpiece</i>	all nebulizers
28140	Nebulizer bulb	all nebulizers

*With elastic

GIMA code	NEBULIZER ACCESSORIES	To be used with
28163	Fasterjet bulb* with Valve System	all nebulizers
28164	Nebjet bulb ** with Double Venturi Effect	all nebulizers
28165	Mouthpiece bi-valve	bulb 28163
28166	Mouthpiece	bulb 28164, 28163
28167	Nasal Prong	bulb 28164, 28163
28168	Facemask-adult***	bulb 28164, 28163
28169	Facemask-pediatry***	bulb 28164, 28163
28145	Hidrofobic, 99%antibacteric filter	28130/1/2/3/4/5

* Standard with Happyneb II, Hospyneb

** Standard with Happyneb I, III, Betaneb

***Without elastic