

STAINLESS STEEL, INSTRUMENT AND DRESSING TROLLEYS



27430

• **27430 INOX TROLLEY**
 AISI 304 Stainless steel trolley with plastic drawer internally divided in 4 sectors. 2 shelves with upstands (3,5 cm) on four sides. Revolving rubber castors Ø 80 mm (two with brake).
 Size: 70x50xh 82 cm.
 Max load (x shelf): 25 kg.
Delivered in kit form.
 Made in Italy.



27433

*Delivered in kit form.
 Made in Italy.*

• **27433 MEDICAZIONE TROLLEY - small**
 Size: 60x40xh 81 cm - weight 8 kg
 • **27434 MEDICAZIONE TROLLEY - medium**
 Size: 70x50xh 82 cm - weight 9 kg
 • **27435 MEDICAZIONE TROLLEY - large**
 Size: 90x60xh 85 cm - weight 15 kg
 AISI 304 stainless steel trolleys, with push handle, 2 shelves with upstands (3,5 cm) on four sides, 4 revolving rubber castors Ø 60/80/90 mm. Max load (x shelf): 25 kg.



27431

• **27431 TER TROLLEY**
 AISI 304 Stainless steel trolley, with push handle. 3 shelves with upstands (3,5 cm) on four sides 4 revolving rubber castors Ø 80 mm.
 Size: 70x50xh 82 cm - weight 13 kg
 Max load (x shelf): 25 kg
Delivered in kit form.
 Made in Italy.



27436

• **27436 DRESSING TROLLEY**
 Dressing trolley made of AISI 304 stainless steel with 2 shelves (70x50 cm), 1 drawer (40x50xh 15 cm), basin holder and basin. 4 revolving castors Ø 100 mm.
 Size: 70x50xh 85 cm - weight: 18 kg.
 Max load (x shelf): 25 kg.
Delivered in kit form. Made in Italy.



45800

• **45800 MEDICAZIONE PLUS TROLLEY**
With removable trays
 Stainless steel frame mounted on castors Ø 80 mm.
 Size: 61x43xh 78 cm. Weight: 10 kg.
 Maximum load (x shelf): 10 kg.
Delivered in kit form.
 Made in Italy.




45805

• **45805 TER PLUS TROLLEY**
With removable trays
 Stainless steel frame, mounted on castors Ø 80 mm.
 Size: 49x30xh 78 cm. Weight: 9 kg.
 Maximum load (x shelf): 10 kg.
Delivered in kit form.
 Made in Italy.

LANCART: STAINLESS STEEL TROLLEY TO BUILD UP

• **45810 LANCART TROLLEY**
 Stainless steel trolley mounted on castors Ø 80 mm. Made in Italy.
 Maximum trolley load: 60 kg.
 You can add following accessories:
 • **45815 SET OF SIDERAILS (for 1 shelf)**
 • **45817 BOTTLE RACK**
 • **45818 WASTE BIN**



Size:
 60x40xh 80 cm

45810



Set of siderails

Bottle rack

Waste bin

45810
 + 2 x 45815
 + 45817
 + 45818

BUILD UP YOUR TROLLEY